


Comparative Peristaltic Pump Tubing Life

The table below depicts hours until failure of 1/4" I.D. x 3/8" O.D. tubing in each corresponding formulation. In each case, a 3-roller pump head was utilized operating at 600 RPM under room temperature (73°F). Tubing failure is measured in hours of use prior to rupture.



The performance of tubing in peristaltic pumping applications is affected by the conditions of use and equipment utilized, along with size and wall thickness of the tubing tested. The data above is presented for information only and should not be utilized for specification purposes.

IMPORTANT: Saint-Gobain Performance Plastics Corporation does not assume any responsibility or liability for any advice furnished by it, for the performance or results of any installation or use of the product or of any final product into which the product may be incorporated by the purchaser and/or user.

It is the user's responsibility to ensure the suitability and safety of Saint-Gobain tubing for all intended uses. Laboratory, field or clinical tests must be conducted in accordance with applicable requirements in order to determine the safety and effectiveness for use of tubing in any particular application.



Saint-Gobain Performance Plastics Corporation
 2664 Gilchrist Road, Akron, OH 44305
 330-798-9240 • 800-798-1554 • 330-798-6968
www.tygon.com