


NEVER MISS A WELL

PlatR is an easy-to-use **tablet application for reliable manual pipetting to microtiter plates and strips**. It is suitable for **96 and 384 well plates**, including white plates and strips.

PlatR lets users import or prepare their own pipetting plans and navigates their pipetting on microtiter plates by **illuminating the wells according to the plan**.

PlatR is an Android application. We are constantly improving its functionalities and offer free software updates.



BENEFITS FOR YOUR LAB

PlatR is a perfect tool for all lab workers from complete beginners to professionals. It can be used in clinical and research laboratories as well as at the universities as a pipetting learning tool for students.



No pipetting mistakes



Full traceability



Increased throughput



Relaxed and precise pipetting



Fast protocol setup



Protocol reusability

[Request a quote](#)


or write us at platr@biosistemika.com for OEM and distribution opportunities.

SAVE WITH PLATR

Task: Preparing plan and pipetting 1 plate per 1 user	Hour/week without PlatR	Hour/week with PlatR	Productivity boost (%)	Weekly savings (\$)	Annual savings (\$)
Preparing pipetting plan	0,5	0,3	83,3	\$13	\$707
Pipetting 96-well plate	0,5	0,4	20,0	\$7	\$353
Pipetting 384-well plate	1,5	1,1	26,7	\$32	\$1.696
Average error rate - 1/96 wells without PlatR	-	-	-	\$9	\$497
Average error rate - 4/384 wells without PlatR	-	-	-	\$4	\$212
96-well plate					\$1.557
384-well plate					\$2.615

*average loaded cost of researcher per hour (80\$) was used for calculations

Annual savings with PlatR (\$)



5 PlatR users that pipette one 96-well plate or one 384-well plate per week can save between \$8,000 and \$13,000 annually.

**ROI
< 6 months**

PLATR PACKAGE

Basic

- > PlatR software unlimited licence
- > Samsung Galaxy Tab
- > Adhesive microtiter plate holders

Optional accessories

- > Stand for adjustable inclination
- > Bluetooth footswitch
- > Holder for PCR plates and strips
- > Holder for ELISA strips
- > Screen protector
- > Touch pen



Request a quote

or write us at platr@biosistemika.com for OEM and distribution opportunities.